

Miscellaneous

Illustration	Description	Notes	Part no.
	Down post, rotatable	<p>For non-height adjustable arms: L 150 mm / 5.9", max. load 12 kg / 26.5 lbs L 300 mm / 11.8", max. load 10 kg / 22 lbs</p> <p>For height adjustable arms: L 150 mm / 5.9", max. load 12 kg / 26.5 lbs L 300 mm / 11.8", max. load 10 kg / 22 lbs</p>	<p>Z-DPR-70.CD_150 Z-DPR-70.CD_300</p> <p>Z-DPR-71.CD_150 Z-DPR-71.CD_300</p>
	Bottle holder	<p>Incl. adapter to mount to CIM wall channel (compatible with GCX channel) Please order CIM wall channel separately</p> <p>Max. bottle diameter: 140 mm / 5.5" Maximum weight: 25 kg / 55.1 lbs (6-fold safety)</p>	Z-GFH-4.1
	Mounting solution for barcode scanner	For Z-bracket right hand side	Z-BCS-Mr.1
	Mounting solution for barcode scanner	Incl. adapter to mount to CIM wall channel (compatible with GCX channel)	Z-BCS-4Z.1